

LTV 1200 Ventilator Presets Table

	Infant (>5 <10 kg)	Pediatric (10-40 kg)	Adult (> 40 kg)
Alarm volume	85	85	85
HP delay	No delay	No delay	No delay
LPP alarm	All breaths	All breaths	All breaths
High rate alarm	80, 30 sec	60, 30 sec	40, 30 sec
High PEEP alarm	5 above set PEEP	5 above set PEEP	5 above set PEEP
Low PEEP alarm	-3 below set PEEP	-3 below set PEEP	-3 below set PEEP
Rise Time	5	4	4
Flow Term	35%	30%	25%
Time Term	0.5	1.0	2.0
PCFT	Off	Off	Off
Leak Comp	On	On	On
Breath rate	20	15	12
Breath type	Pressure	Pressure	Volume
Volume	50	250	500
Inspiratory Time	0.3	0.7	1.0
Pressure Control	15	15	15
Pressure Support	10	10	10
Sensitivity	2	3	3
High Pressure	30	30	40
Low Pressure	10	10	10
Low Min Volume	0.5 ml	1.0 ml	3.0 ml
PEEP	0	0	0
Mode	Assist/Control	Assist/Control	Assist/Control